

Experiment-1

Name: M Jaswanth

Regno: 20bcd7171

Title: Determination of molarity of HCl acid

Aim: To determine the molarity of the given dil. HCl acid by P-H metry by using 0.1M NaOH solution.

Procedure: The standard NaOH solution is taken in 25ml burette (well cleaned). Then 20ml of dil HCl solution will be taken in 250ml of beaker and then mixed with 80ml of deionised water. Then PH-electrode is placed in beaker to measure the PH of the solution. Now by adding NaOH drop to drop it means by starting with 1ml is added and P-H is measured such that NaOH is added to the solution with addition of 1ml by successive additions up to 15ml.


Observation: P-H will increase initially from range of acid at $\text{PH} = 2-4$ and it sharply increases at near $\text{PH} = 7$. Not more than 7. If we add ~~PH~~^{NaOH} like that only then after $\text{PH} = 7$ the PH will be increased to $\text{PH} = 9-11$ that is alkaline range.

S NO	Addition of NaOH (cm ³)	Volume of NaOH (cm ³)	p ^H measured
1	0	0	1.74
2	0.5	0.5	1.75
3	0.5	1	1.78
4	0.5	1.5	1.8
5	0.5	2	1.82
6	0.5	2.5	1.85
7	0.5	3	1.88
8	0.5	3.5	1.91
9	0.5	4	1.94
10	0.5	4.5	1.99
11	0.5	5	2.04
12	0.5	5.5	2.09
13	0.5	6	2.12
14	0.5	6.5	2.15
15	0.5	7	2.17
16	0.5	7.5	2.31
17	0.5	8	2.42
18	0.5	8.5	2.59
19	0.5	9	2.84
20	0.5	9.5	3.84
21	0.5	10	9.48
22	0.5	10.5	10.38
23	0.5	11	10.65
24	0.5	11.5	10.78
25	0.5	12	10.88
26	0.5	12.5	10.95
27	0.5	13	11.01
28	0.5	13.5	11.06
29	0.5	14	11.09

Scale:

X-axis 1cm = 1unit

Y-axis 1cm = 1unit



Volume of NaOH

Calculations:

V_1 = Volume of HCl given = 20ml.

M_1 = Strength of HCl solution = M

V_2 = volume of NaOH solution at end point = 9.8ml

M_2 = Strength of NaOH solution = 0.1 M.

$$\text{Strength of HCl solution } M_1 = \frac{V_2 \times M_2}{V_1}$$

$$= \frac{9.8 \times 0.1}{20}$$

$$= 0.049 \text{ M}$$

$$\text{Amount of HCl in 1L of given solution} = 0.049 \times 36.45$$

$$= 1.786 \text{ g/L}$$

Results:

Strength of given HCl solution = 0.049 M

Amount of HCl present in 1L = 1.786 g/L

Answers for the questions:

1) A) Solid NaOH consists of Na^+ and OH^- ions packed into crystalline lattice. When solid is added to water the ions float leading apart to extra OH^- ions in the water



The resulting large concentration of OH^- makes the solution more basic and leads to a dramatic increase in pH.

2) A) When an acidic solution is diluted with water the concentration of H^+ ions decreases and pH of the solution increases towards 7. To make the pH change by 1, a tenfold dilution is required. The acid is becoming less acidic

3) A) Acetic acid is a weak electrolyte. It is not completely ionized and hence gives less H^+ ion concentration, HCl is a strong acid. It is completely ionized giving more H^+ ion concentration. As $\text{pH} = -\log[\text{H}^+]$, less the $[\text{H}^+]$ greater will be pH