Estimation of water hardness

Volume of Hard water taken = 10 ml

S.No	Vol. of hard	Burette reading		Volume of 0.01 M	Average volume of
	water	Initial	Final	EDTA (ml)	EDTA (ml)
	sample (ml)				
1	10	0	11.5	11.5	
2	10	0	11.7	11.7	
3	10	0	11.4	11.4	

Questions: (Please answer the following questions in 2 or 3 sentences)

- 1. Why does colour of solution changes from wine red to blue at end point?
- 2. What is the role of buffer in this experiment?
- 3. Why is EDTA used in this titration?